





**TERUMO® Surshield™ Safety Winged Infusion Set**

**<INSTRUCTIONS FOR USE: IV ADMINISTRATION>**

Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on ALL patients.

**PRECAUTION:** Keep hands behind the needle at all times during use and disposal.

**SET UP**

1. Open the package.
2. Connect an administration set, syringe or etc. and secure.
3. Proceed with priming and confirm that administration fluid comes from the needle.

**VENIPUNCTURE AND ADMINISTRATION**

4. Flip the safety shield back away from the needle toward the tubing. Grasp wings.
5. Remove the needle protector. Caution: Care should be taken not to touch the needle.
6. Perform venipuncture and confirm proper positioning in the vessel.
7. Secure the position of the winged infusion set per facility protocol.

**AFTER USE**

8. Remove tape, if present, from wings.
9. Flip the safety shield forward toward needle. Grasp the wing and shield between your thumb and finger. Completely remove the needle from the puncture site and apply digital pressure to the site using a sterile gauze pad. (Fig. 1)
10. While maintaining the same grip on the device (per step 9 above), activate the safety feature by performing either of the following methods:  
 Recommended (one-handed) method: (Reference OSHA 1910.1030 (d)(2)(vii)(B))  
 Pinch your finger and thumb together (Fig. 2) (or press the safety shield against a hard surface such as a bedside table) until an audible click is heard and the safety feature is visually confirmed to be activated. (Fig. 3)  
**OR**  
 Alternate (two-handed) method:  
 While keeping both hands behind the needle, grasp one of the wing tabs with your free hand to stabilize the device. (Fig. 4) (If necessary, reposition your original grip per step 9 being careful not to press against an exposed area of the needle.) Then, pinch your finger and thumb together (Fig. 2) until an audible click is heard and the safety feature is visually confirmed to be activated. (Fig. 3) Caution: Needle must be carefully and firmly placed in the lock position of the safety shield.  
 11. Do not attempt to deactivate the safety device by separating the needle from the safety shield.  
 12. Dispose of used needles and materials following the policies and procedures of your facility, as well as federal and local regulations for "Sharps Disposal."

**<INSTRUCTIONS FOR USE: BLOOD COLLECTION>**

Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on ALL patients.

**PRECAUTION:** Keep hands behind the needle at all times during use and disposal.

**SET UP**

1. Open the package.
2. Connect a syringe, luer adapter, or other blood collection device and secure.

**BLOOD COLLECTION**

3. Flip the safety shield back away from the needle toward the tubing. Grasp wings.
4. Remove the needle protector. Caution: Care should be taken not to touch the needle.
5. Perform venipuncture following facility protocol.

**AFTER USE**

6. Remove tape, if present, from wings.
7. Flip the safety shield forward toward needle. Grasp the wing and shield between your thumb and finger. Completely remove the needle from the puncture site and apply digital pressure to the site using a sterile gauze pad. (Fig. 1)
8. While maintaining the same grip on the device (per step 7 above), activate the safety feature by performing either of the following methods:  
 Recommended (one-handed) method: (Reference OSHA 1910.1030 (d)(2)(vii)(B))  
 Pinch your finger and thumb together (Fig. 2) (or press the safety shield against a hard surface such as a bedside table) until an audible click is heard and the safety feature is visually confirmed to be activated. (Fig. 3)  
**OR**  
 Alternate (two-handed) method:  
 While keeping both hands behind the needle, grasp one of the wing tabs with your free hand to stabilize the device. (Fig. 4) (If necessary, reposition your original grip per step 9 being careful not to press against an exposed area of the needle.) Then, pinch your finger and thumb together (Fig. 2) until an audible click is heard and the safety feature is visually confirmed to be activated. (Fig. 3) Caution: Needle must be carefully and firmly placed in the lock position of the safety shield.  
 9. Do not attempt to deactivate the safety device by separating the needle from the safety shield.  
 10. Dispose of used needles and materials following policies and procedures of your facility, as well as federal and local regulations for "Sharps Disposal."

®: Registered Trademark Surshield is a trademark of TERUMO CORPORATION.

Rev. Date: 08/01