

# Introducing Sur-Vet<sup>®</sup> Secure



Topical adhesive for catheter  
securement and site protection



**Secure catheters. Fewer complications.<sup>1</sup>**



**TERUMO**

For additional information, visit [TerumoTMP.com](http://TerumoTMP.com)



**Sur-Vet<sup>®</sup> Secure**  
Catheter Securement Adhesive

# A new approach to vascular access device securement.

## Simple, strong and safe.

Sur-Vet® Secure provides benefits and advantages over traditional methods.

- **Up to 3X stronger** securement than transparent film dressings alone and superior adhesiveness compared with mechanical securement devices<sup>2</sup>
  - Reduces catheter movement, migration, and dislodgement<sup>1</sup>
- **Effectively seals** the catheter insertion site
  - Can help reduce early dressing changes and eliminate 24-48 hour initial dressing change protocols<sup>3</sup>
  - Formulation reported to exhibit activity against gram-negative and gram-positive bacteria, yeast, and fungi eliminating greater than 8-logs after 3 minutes of contact in in-vitro studies<sup>4</sup>
  - Allows the direct visualization of insertion site at all times, unlike some antimicrobial dressings

### Sur-Vet® Secure Catheter Securement Adhesive

Description	Product Code	Fill Volume	Packaging
Sur-Vet Secure Catheter Securement Adhesive	SV-300V	2.5ml	50 applications per bottle



Sur-Vet Secure adhesive can be used on peripheral IVs by placing 2-3 drops over the insertion site and underneath the catheter hub.

## FAQ's:

- **How many drops are required to secure a catheter? Where should the drops be placed?**
  - Apply one drop of adhesive directly on top of the catheter injection site. An additional 1-2 drops of adhesive can be used under the catheter hubs and/or wings for a stronger securement.
- **How is the Sur-Vet Secure adhesive removed?**
  - Standard adhesive remover or acetone is recommended. If removal is not required early, Sur-Vet Secure will naturally wear away, typically within 5-7 days, for easy removal of catheter.
- **How many catheter placements can one bottle yield?**
  - Sur-Vet Secure can be used to place approximately 50 Sur-Vet SurFlo catheters.
- **Will Sur-Vet Secure work to seal wounds?**
  - There are no claims to support wound closure at this time.

©2018 Adhezion Biomedical 08/18. All rights reserved. SPI-SS01-1801

1. Intravascular canine model 2017, reduced dislodgement, data on file.

2. Data on file. Average securement strength calculated from measured values. Nov/Dec 2016.

3. Kleidon, et al. A Pilot Randomized Controlled Trial of Novel Dressing and Securement Techniques in 101 Pediatric Patients. J Vasc Interv Radiol. 2017 Sep 18.

4. Prince, et al. Immobilization and Death of Bacteria by Flora Seal® Microbial Sealant. International Journal of Pharmaceutical Science Invention, 2017.

**RX ONLY.** The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

TERUMO, TERUMO, SUR-VET and SURFLO are trademarks owned by Terumo Corporation, Tokyo, Japan, and are registered with the U.S. Patent and Trademark Office. ©2020 Terumo Medical Corporation 07/20. All rights reserved. PM-03416

TerumoTMP.com



TERUMO MEDICAL PRODUCTS