

TERUMO® Conventional Needles & Syringes




Superior sharpness and smoothness time after time

- Ultra Sharp Needles
- Clear Barrels
- Bold, Precise Scale Markings
- Smooth Plunger Action
- Tamper Evident Peel Packs

When comfort, performance and cost effectiveness are important, choose Terumo. Our ultra-thin wall needles are precisely honed and beveled for ease of penetration and less patient trauma. Clear syringe barrels with bold scale markings permit quick visual inspection for precise dosage. An integrated plunger stop helps prevent the plunger from being accidentally pulled out. This product contains no components made of natural rubber latex. And our tamper evident blister packs result in minimal waste disposal.

All Competitors



* In Independent laboratory testing, Terumo needles were shown, on average, to be the sharpest among brands tested:

10% sharper than BD (22G - 30G)
16% sharper than Covidien (22G - 27G)

When you think hypodermics...
Think sharp. Think Terumo.



Terumo® Needles & Syringes

TERUMO Hypodermic Needles

PACKED: 100/box, 10 boxes/case (1,000)

Product Code	Description
NN-1825R	18G x 1", Thin Wall Needle
NN-1838R	18G x 1 1/2", Thin Wall Needle
NN-1925R	19G x 1", Thin Wall Needle
NN-1938R	19G x 1 1/2", Thin Wall Needle
NN-2025R	20G x 1", Ultra-Thin Wall Needle
NN-2038R	20G x 1 1/2", Ultra-Thin Wall Needle
NN-2125R	21G x 1", Ultra-Thin Wall Needle
NN-2138R	21G x 1 1/2", Ultra-Thin Wall Needle
NN-2225R	22G x 1", Ultra-Thin Wall Needle
NN-2238R	22G x 1 1/2", Ultra-Thin Wall Needle
NN-2325R	23G x 1", Ultra Thin Wall Needle
NN-2338R	23G x 1 1/2", Ultra-Thin Wall Needle
NN-2516R	25G x 5/8", Ultra Thin Wall Needle
NN-2522R	25G x 7/8", Ultra Thin Wall Needle
NN-2525R	25G x 1", Ultra Thin Wall Needle
NN-2538R	25G x 1 1/2", Regular Wall Needle
NN-2613R	26G x 1/2", Regular Wall Needle
NN-2713R	27G x 1/2", Regular Wall Needle
NN-2732R	27G x 1 1/4", Regular Wall Needle
NN-3013R	30G x 1/2", Regular Wall Needle
NN-3025R	30G x 1", Regular Wall Needle

TERUMO Tuberculin Syringes

PACKED: 100/box, 10 boxes/case (1,000)

Product Code	Description
SS-01T	1cc TB Syringe without needle, Graduations: 0.01cc (1cc total)
SS-01T2516	1cc TB Syringe with 25G x 5/8" removable needle, Graduations: 0.01 cc (1cc total)
SS-01T2609	1cc TB Syringe with 26G x 3/8" removable needle, Graduations: 0.01 cc (1cc total)
SS-01T2713	1cc TB Syringe with 27G x 1/2" removable needle, Graduations: 0.01cc (1cc total)

SurSaver® 1CC Syringe with Low Dead Space

PACKED: 100/box, 5 boxes/case (500)

Product Code	Description
SS01D2313	1cc syringe with 23G x 1/2" fixed needle
SS01D2516	1cc syringe with 25G x 5/8" fixed needle
SS01D2713	1cc syringe with 27G x 1/2" fixed needle

TerumoTMP.com/Hypodermics



and TERUMO are registered trademarks owned by Terumo Corporation, Tokyo, Japan. SurSaver is a registered trademark owned by Terumo Medical Corporation, Elkton, MD, U.S.A. TMP-0158-11192015



TERUMO Hypodermic Syringes Without Needle

Product Code Description

Packed: 100/box, 10 boxes/case (1,000)

SS-03L	3cc Luer Lock Tip Syringe without Needle, Graduations: 0.1cc (3cc total)
SS-03S	3cc Luer Slip Tip Syringe without Needle, Graduations: 0.1cc (3cc total)
Packed: 100/box, 6 boxes/case (600)	
SS-05L	5cc Luer Lock Tip Syringe without Needle, Graduations: 0.2cc (6cc total)
SS-05S	5cc Luer Slip Tip Syringe without Needle, Graduations: 0.2cc (6 cc total)
SS-10L	10cc Luer Lock Tip Syringe without Needle, Graduations: 0.2cc (12cc total)
SS-10S	10cc Slip Tip Syringe without Needle, Graduations: 0.2cc (12cc total)
SS-10ES	10cc Eccentric Luer Tip Syringe without Needle, Graduations: 0.2cc (12cc total)
Packed: 50/box, 4 boxes/case (200)	
SS-20L2	20cc Luer Lock Tip Syringe without Needle, Graduations: 1.0cc (25cc total)
SS-20S2	20cc Luer Slip Tip Syringe without Needle, Graduations: 1.0cc (25cc total)
SS-20ES2	20cc Eccentric Luer Tip Syringe without Needle, Graduations: 1.0cc (25cc total)
Packed: 25/box, 8 boxes/case (200)	
SS-30L	30cc Luer Lock Tip Syringe without Needle, Graduations: 1.0cc (35cc total)
SS-30S	30cc Luer Slip Tip Syringe without Needle, Graduations: 1.0cc (35cc total)
SS-30ES	30cc Eccentric Luer Tip Syringe without Needle, Graduations: 1.0cc (35cc total)
Packed: 25/box, 4 boxes/case (100)	
SS-60L	60cc Luer Lock Tip Syringe without Needle, Graduations: 1.0cc (60cc total)
SS-60ES	60cc Eccentric Luer Tip Syringe without Needle, Graduations: 1.0cc (60cc total)
SS-60C	60cc 2oz. Catheter Tip Syringe without Needle, Graduations: 1.0cc (60cc total)

TERUMO Hypodermic Syringes With Needle

Product Code Description (Luer Lock Tip Syringe)

Packed: 100/box, 10 boxes/case (1,000)

SS-03L2025	3cc with 20G x 1" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2038	3cc with 20G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2125	3cc with 21G x 1" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2138	3cc with 21G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2225	3cc with 22G x 1" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2238	3cc with 22G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2325	3cc with 23G x 1" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2338	3cc with 23G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2516	3cc with 25G x 5/8" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
SS-03L2525	3cc with 25G x 1" Ultra Thin Wall Needle, Graduations: 0.1cc (3cc total)
Packed: 100/box, 6 boxes/case (600)	
SS-05L2025	5cc with 20G x 1" Ultra Thin Wall Needle, Graduations: 0.2cc (6cc total)
SS-05L2038	5cc with 20G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.2cc (6cc total)
SS-05L2125	5cc with 21G x 1" Ultra Thin Wall Needle, Graduations: 0.2cc (6cc total)
SS-05L2138	5cc with 21G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.2cc (6cc total)
SS-05L2225	5cc with 22G x 1" Ultra Thin Wall Needle, Graduations: 0.2cc (6cc total)
SS-05L2238	5cc with 22G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.2cc (6cc total)
SS-10L2025	10cc with 20G x 1" Ultra Thin Wall Needle, Graduations: 0.2cc (12cc total)
SS-10L2038	10cc with 20G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.2cc (12cc total)
SS-10L2125	10cc with 21G x 1" Ultra Thin Wall Needle, Graduations: 0.2cc (12cc total)
SS-10L2138	10cc with 21G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.2cc (12cc total)
SS-10L2238	10cc with 22G x 1 1/2" Ultra Thin Wall Needle, Graduations: 0.2cc (12cc total)

Pharmacy Trays

Product Code Description (Luer Lock Tip without Needle)

SS10LT	10cc, Graduations: 0.2cc (12cc total) Packed: 25 syringes/tray, 8 trays/case (200)
SS20LT	20cc, Graduations: 1.0cc (25cc total) Packed: 25 syringes/tray, 8 trays/case (200)
SS30LT	30cc, Graduations: 1.0cc (35cc total) Packed: 10 syringes/tray, 12 trays/case (120)
SS60LT	60cc, Graduations: 1.0cc (60cc total) Packed: 10 syringes/tray, 12 trays/case (120)

TERUMO MEDICAL PRODUCTS