

TERUMO

Product Catalog (Canada)

Conventional Needles · I.V. Catheters · Micro-Collection Products · Safety Products · Winged Infusion Sets



Table of Contents

Hypodermics

SurGuard®3 Safety Hypodermic Needle	3-4
Terumo Standard Hypodermic Needles	5
Terumo Tuberculin Syringes	6
Terumo Hypodermic Syringes with Needle	6
Terumo Hypodermic Syringes without Needle	6

Micro-collection

Surshield® Safety Blood Collection Sets	7
Multi-Sample Luer Adapter	7
Capiject® Safety Lancets	8

Winged Infusion Needle Sets


Surshield® Safety Winged Infusion Sets	9
SURFLO® Winged Infusion Sets	9

IV Catheters

SURFLO® ETFE I.V. Catheters	10
SurFlash® Polyurethane I.V. Catheters	10

Frequently Asked Questions and Technical Support

11



Pick how you click.




With 3 ways to activate and a broader range of sizes than most other brands, SurGuard®3 is the industry's most adaptable safety needle. Add separate thumb and finger grips, sharper needles, and competitive pricing, and you'll understand why we say SurGuard®3 is the one needle that fits *every* clinician's style.

- 3 modes of activation—thumb, finger, or hard surface
- Separate finger and thumb grips with ridges and nonslip surfaces
- Hub and needle locks for security post activation
- Broad range of sizes
- Firm injection feel
- Sharpest needles among major brands tested*
 - 10% sharper than similar needles from Becton, Dickinson and Company
 - More comfortable injections for patients

Bottom line: No matter what you're looking for in a needle, you're looking for SurGuard®3.

*Reference: Data on file. Terumo Medical Products, April 2009.

Order the industry's most adaptable safety needle today!



SurGuard®3 Safety Hypodermic Needle

Packaging: 100/box, 8 boxes/case (800)

Order No.	Description
SG3-1825	18G x 1" needle
SG3-1838	18G x 1 1/2" needle
SG3-2025	20G x 1" needle
SG3-2038	20G x 1 1/2" needle
SG3-2051	20G x 2" needle
SG3-2125	21G x 1" needle
SG3-2138	21G x 1 1/2" needle
SG3-2225	22G x 1" needle
SG3-2238	22G x 1 1/2" needle
SG3-2325	23G x 1" needle
SG3-2338	23G x 1 1/2" needle
SG3-2516	25G x 5/8" needle
SG3-2525	25G x 1" needle
SG3-2538	25G x 1 1/2" needle
SG3-2713	27G x 1/2" needle
SG3-3013	30G x 1/2" needle

SurGuard®3 Hypodermic Syringes with Safety Needle

Packaging: 100/box, 4 boxes/case (400)

Order No.	Description
SG3-01T2516	1CC Syringe with 25G x 5/8" needle
SG3-01T2713	1CC Syringe with 27G x 1/2" needle
SG3-03L2025	3CC Syringe with 20G x 1" needle
SG3-03L2038	3CC Syringe with 20G x 1 1/2" needle
SG3-03L2125	3CC Syringe with 21G x 1" needle
SG3-03L2138	3CC Syringe with 21G x 1 1/2" needle
SG3-03L2225	3CC Syringe with 22G x 1" needle
SG3-03L2238	3CC Syringe with 22G x 1 1/2" needle
SG3-03L2325	3CC Syringe with 23G x 1" needle
SG3-03L2516	3CC Syringe with 25G x 5/8" needle
SG3-03L2525	3CC Syringe with 25G x 1" needle



You know what you like.
And *whatever* your preference,
you'll find it in SurGuard®3!



SurGuard®3
SAFETY HYPODERMIC NEEDLE

Our double-bevel design is precision-ground and honed to exacting quality standards. This attention to detail enables Terumo needles to penetrate skin with less force, which can mean greater comfort for patients.




Terumo® Standard Hypodermic Needles

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
■ NN1825R	18G x 1" TW Needle
■ NN1838R	18G x 1½" TW Needle
■ NN1925R	19G x 1" TW Needle
■ NN1938R	19G x 1½" TW Needle
■ NN2025R	20G x 1" UTW Needle
■ NN2038R	20G x 1½" UTW Needle
■ NN2125R	21G x 1" UTW Needle
■ NN2138R	21G x 1½" UTW Needle
■ NN2225R	22G x 1" UTW Needle
■ NN2238R	22G x 1½" UTW Needle
■ NN2325R	23G x 1" UTW Needle
■ NN2338R	23G x 1½" UTW Needle
■ NN2516R	25G x 5/8" UTW Needle
■ NN2525R	25G x 1" UTW Needle
■ NN2538R	25G x 1½" RW Needle
■ NN2613R	26G x ½" RW Needle
■ NN2713R	27G x ½" RW Needle
■ NN2732R	27G x ¼" RW Needle
■ NN3013R	30G x ½" RW Needle
■ NN3025R	30G x 1" RW Needle

RW - Regular Wall
 TW - Thin Wall
 UTW - Ultra Thin Wall

All Competitors



* In Independent laboratory testing, Terumo needles were shown, on average, to be the sharpest among brands tested:

10% sharper than BD (22G - 30G)

16% sharper than Covidien (22G - 27G)

Conventional Syringes

1cc Tuberculin Syringes

Graduations: 0.01cc (1cc total)

TERUMO TB Syringes have clear, bold scale markings which assure dosage control and accuracy. The package is tamper evident, for greater protection of the sterile contents.

Packed: 100/box, 10 boxes/case (1,000)

Product Code Description

SS-01T	1cc TB Syringe without needle, with luer slip
SS-01T2516	1cc TB Syringe with 25G x 5/8" removable needle
SS-01T2609	1cc TB Syringe with 26G x 3/8" removable needle
SS-01T2713	1cc TB Syringe with 27G x 1/2" removable needle



Catheter Tip



Luer Lock Tip



Eccentric Luer Tip



Luer Slip Tip

3cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
SS-03L2025	3cc, 20G x 1" UTW
SS-03L2038	3cc, 20G x 1 1/2" UTW
SS-03L2125	3cc, 21G x 1" UTW
SS-03L2138	3cc, 21G x 1 1/2" UTW
SS-03L2225	3cc, 22G x 1" UTW
SS-03L2238	3cc, 22G x 1 1/2" UTW
SS-03L2325	3cc, 23G x 1" UTW
SS-03L2338	3cc, 23G x 1 1/2" UTW
SS-03L2516	3cc, 25G x 5/8" UTW
SS-03L2525	3cc, 25G x 1" UTW

5cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.2cc (6cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code	Description
SS-05L2025	5cc, 20G x 1" UTW
SS-05L2038	5cc, 20G x 1 1/2" UTW
SS-05L2125	5cc, 21G x 1" UTW
SS-05L2138	5cc, 21G x 1 1/2" UTW
SS-05L2225	5cc, 22G x 1" UTW
SS-05L2238	5cc, 22G x 1 1/2" UTW

10cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.2cc (12cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code	Description
SS-10L2025	10cc, 20G x 1" UTW
SS-10L2038	10cc, 20G x 1 1/2" UTW
SS-10L2125	10cc, 21G x 1" UTW
SS-10L2138	10cc, 21G x 1 1/2" UTW

UTW - Ultra Thin Wall

3cc Syringes without Needle

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

Product Code Description

SS-03L	Luer Lock Tip
SS-03S	Luer Slip Tip

5cc Syringes without Needle

Graduations: 0.2cc (6cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code

Product Code	Description
SS-05L	5cc Luer Lock Tip
SS-05S	5cc Luer Slip Tip

10cc Syringes without Needle

Graduations: 0.2cc (12cc total)

Packed: 100/box, 6 boxes/case (600)

Product Code

Product Code	Description
SS-10L	10cc Luer Lock Tip
SS-10S	10cc Slip Tip
SS-10ES	10cc Eccentric Luer Tip

20cc Syringes without Needle

Graduations: 1.0cc (25cc total)

Packed: 50/box, 4 boxes/case (200)

Product Code

Product Code	Description
SS-20L2	20cc Luer Lock Tip
SS-20S2	20cc Luer Slip Tip
SS-20ES2	20cc Eccentric Luer Tip

30cc Syringes without Needle

Graduations: 1.0cc (35cc total)

Packed: 25/box, 8 boxes/case (200)

Product Code

Product Code	Description
SS-30L	30cc Luer Lock Tip
SS-30S	30cc Luer Slip Tip
SS-30ES	30cc Eccentric Luer Tip

60cc Syringes without Needle

Graduations: 1.0cc (60cc total)

Packed: 25/box, 4 boxes/case (100)

Product Code

Product Code	Description
SS-60L	60cc Luer Lock Tip
SS-60ES	60cc Eccentric Luer Tip
SS-60C	60cc 2oz. Catheter Tip

These products have no components made of natural rubber latex.



Surshield® Safety Winged Blood Collection Sets

Packed: 50/box, 10 boxes/case (500)

Product Code Description

(12" Tubing)

MN*SVS21B30	21G x 3/4" needle
MN*SVS23B30	23G x 3/4" needle
MN*SVS25B30	25G x 3/4" needle

(7" Tubing)

MN*SVS21B18	21G x 3/4" needle
MN*SVS23B18	23G x 3/4" needle



Multi-Sample Luer Adapter

Sterile and disposable, adapters fit into standard infusion sets transforming them into winged blood collection sets. This allows multiple draws into evacuated tubes. Accepts standard tube holders.

Packed: 100/box, 30 boxes/case (3,000)

Product Code Description

MN*2000TE	VENOJECT® Multi-Sample Luer Adapter (male).
-----------	---




Blood Collection Sets

Safety


Terumo's Surshield Safety Winged Needle Sets with one-handed activation provides a solution with a cost-effective, easy to use design. Press shield over the needle until it locks into place.

- Unique needle locking mechanism helps reduce the risk of needlestick injuries
- One-handed activation for safety and ease of use
- No change in needle insertion technique
- Simple, built-in safety design is easy to use

1. Flip safety shield toward needle. Remove the needle and apply digital pressure to the site.



2. With the wing and shield between your thumb and index finger, press the safety shield against a hard surface such as a bedside table until an audible click is heard.



3. Visually confirm activation of the safety feature. (10/08)



Micro-Collection

Capiject Safety Lancets

High quality, single-use lancets offer the form, fit and function you've come to rely on

- **Never see or handle the sharp:**
Automatic blade or needle retraction styles - the sharp edge is enclosed before and after use
- **No special staging required:**
Ready to use right out of the box
- **One-step activation:**
Convenient activation for precise and consistent results
- **Color-coded sizes:**
Quick identification of the right product

Packed: 200/box, 10 boxes/case (2,000)



Catalog Number	Color	Description	Width/Gauge	Depth
200101	Pink	Blade	1.5mm	1.0mm
200102	Green	Blade	1.5mm	1.5mm
200103	Blue	Blade	1.5mm	2.0mm
200104	Orange	Needle	23G	2.25mm
200105	Purple	Needle	28G	1.25mm

These products have no components made of natural rubber latex.



**Surshield® Safety
Winged Infusion Sets**

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
(3½" Tubing)	
■ SV*S19BLS	19G x ¾" Needle
■ SV*S21BLS	21G x ¾" Needle
■ SV*S23BLS	23G x ¾" Needle
■ SV*S25BLS	25G x ¾" Needle

(12" Tubing)	
■ SV*S19BL	19G x ¾" Needle
■ SV*S21BL	21G x ¾" Needle
■ SV*S23BL	23G x ¾" Needle
■ SV*S25BL	25G x ¾" Needle



SURFLO® Winged Infusion Sets

Product Code	Description
(3½" Tubing)	
Packed: 100/box, 5 boxes/case (500)	
■ SV*19BLS	19G x ¾" TW
■ SV*21BLS	21G x ¾" UTW
■ SV*23BLS	23G x ¾" UTW
■ SV*25BLS	25G x ¾" TW

(8" Tubing)	
Packed: 100/box, 10 boxes/case (1,000)	
■ SV*25EL	25G x ½" TW
■ SV*27EL	27G x ½" RW

(12" Tubing)	
Packed: 50/box, 10 boxes/case (500)	
■ SV-18BLK	18G x ¾" TW
■ SV-19BLK	19G x ¾" TW
■ SV-21BLK	21G x ¾" UTW
■ SV-22BLK	22G x ¾" UTW
■ SV-23BLK	23G x ¾" UTW
■ SV-25BLK	25G x ¾" TW

RW - Regular Wall
TW - Thin Wall
UTW - Ultra Thin Wall

I.V. Catheters



Terumo's Polyurethane I.V. Catheter enables quick and easy confirmation of placement in the vein.

With SurFlash, visual confirmation of catheter placement occurs due to a unique flashback feature. Flashback appears between the specially grooved needle and the catheter, providing medical staff with a higher certainty of proper catheter placement. The specially grooved needle also improves the ability to visualize flashback in patients where venipuncture is more difficult, and because SurFlash is made of polyurethane, which is firm at room temperature and softens at body temperature, patient comfort is enhanced, without compromising proper technique.

SURFLO® ETFE I.V. Catheters (Large Photo)

A sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue. Preview chamber in clear sure-grip hub gives immediate indication of proper placement. Flexible thinwall catheter with large inner diameter assures good blood flow. All this plus safe, easy-to-handle, tamper-evident packaging make SURFLO® I.V. Catheters a great clinical value.

Packed: 50/box, 4 boxes/case (200)

Product Code	Description
SR-OX1451CA	14G x 2"
SR-OX1632CA	16G x 1 1/4"
SR-OX1651CA	16G x 2"
SR-OX1664CA	16G x 2 1/2"
SR-OX1832CA	18G x 1 1/4"
SR-OX1851CA	18G x 2"
SR-OX1864CA	18G x 2 1/2"
SR-OX2025CA	20G x 1"
SR-OX2032CA	20G x 1 1/4"
SR-OX2051CA	20G x 2"
SR-OX2225CA	22G x 1"
SR-OX2419CA	24G x 3/4"

SurFlash® Polyurethane I.V. Catheters (Not Pictured)

- Visual confirmation of catheter placement during venipuncture occurs in a flash!
- Made of polyurethane, which makes it extremely flexible, kink-resistant and durable
- Provides medical staff with higher certainty of catheter placement due to unique flashback feature
- Improves visualization for difficult venipunctures

Packed: 50/box, 4 boxes/case (200)

Product Code	Description
SR*FF1451	14G x 2"
SR*FF1651	16G x 2"
SR*FF1664	16G x 2 1/2"
SR*FF1832	18G x 1 1/4"
SR*FF1851	18G x 2"
SR*FF1864	18G x 2 1/2"
SR*FF2025	20G x 1"
SR*FF2032	20G x 1 1/4"
SR*FF2051	20G x 2"
SR*FF2225	22G x 1"
SR*FF2419	24G x 3/4"

These products have no components made of natural rubber latex.

Frequently Asked Questions and Technical Support

Who is the Sales Rep for my area?

Just log-on to www.TerumoTMP.com and select “Find Your Rep” from our homepage.

How can I obtain samples?

Need samples? We’ve now made it faster and easier for you to obtain samples of our products through your local TMP rep. Just log-on to www.TerumoTMP.com and select “Find Your Rep” from our homepage.

How can I order Terumo products through a distributor?

Please contact your authorized distributor. To find a distributor, visit us at www.TerumoTMP.com/Support to find a List of Authorized Terumo Distributors. Distributors should call Customer Administration for all questions related to orders and invoices: 800.888.3786

How can I obtain pricing?

You must check with your local distributor for pricing since our business is conducted through authorized distributors.

What is the shelf life of Terumo Products?

Standard & Safety Syringes & Needles: 5 years
Standard IV Catheters: 3 years
Standard & Safety Blood Collection/Infusion: 3 years
Capject Lancets: 5 years

What’s the difference between SurFlash® & Surflo® IV Catheters?

SurFlash is an I.V. Polyurethane Catheter which allows for a visual confirmation of catheter placement due to its unique flashback feature. SURFLO is our regular ETFE I.V. Catheter. Both are radio-opaque and these products have no components made of natural rubber latex.

Where can I get a Certificate of Analysis on products I have ordered?

Visit us at <http://www.terumomedical.com/Default.asp?dest=Certification/MainCertMenu.asp>, choose your product category and enter your lot number and receive a Certificate of Analysis.

What does the “Think green. Think Terumo.” Icon represent?

In 1999, Terumo Corporation developed a Basic Environmental Policy as part of our corporate mission of “Contributing to Society through Healthcare.” Since then, Terumo has engaged in a range of activities aimed at protecting the global environment. Terumo’s green products are Terumo Syringes, Surshield Safety Winged Blood Collection Sets, Surshield Safety Winged Infusion Sets, Surflo Blood Collection Sets and all Survet Hypodermic Syringes.

www.TerumoTMP.com

 **TERUMO** and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. CAPIJECT and SURGUARD are registered trademarks of Terumo Medical Corporation, Somerset, NJ. SURFLO, SURSHIELD and VENOJECT are registered trademarks of Terumo Corporation, Tokyo, Japan and they are registered with the U.S. Patent and Trademark Office. Approval # TMP-0319-10112017