

SURFLO® I.V. Catheters



SURFLO® ETFE I.V. Catheters

A sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue. Preview chamber in clear sure-grip hub gives immediate indication of proper placement. Flexible thinwall catheter with large inner diameter assures good blood flow. All this plus easy-to-handle, tamper-evident packaging make SURFLO® I.V. Catheters a great clinical value.

Packed: 50/box, 4 boxes/case (200)

Product Code Description

■ SR-OX1451CA	14G x 2"
■ SR-OX1632CA	16G x 1¼"
■ SR-OX1651CA	16G x 2"
■ SR-OX1664CA	16G x 2½"
■ SR-OX1832CA	18G x 1¼"
■ SR-OX1851CA	18G x 2"
■ SR-OX1864CA	18G x 2½"
■ SR-OX2025CA	20G x 1"
■ SR-OX2032CA	20G x 1¼"
■ SR-OX2051CA	20G x 2"
■ SR-OX2225CA	22G x 1"
■ SR-OX2419CA	24G x ¾"



Call your Terumo representative at **1.800.888.3786**
or visit www.TerumoTMP.com and click on Find Your Local Terumo Rep.

TerumoTMP.com

These products contain no components made of natural rubber latex.

RX ONLY. The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURFLO is a trademark owned by Terumo Corporation, Tokyo, Japan. PM-03446

