



TERUMO

Product Catalog (Canada)

Conventional Needles · I.V. Catheters · Blood Collection Sets · Safety Products · Winged Infusion Sets



Table of Contents

Hypodermics

| | |
|---|-----|
| SurGuard®3 Safety Hypodermic Needle | 3-4 |
| Terumo Standard Hypodermic Needles | 5 |
| Terumo Tuberculin Syringes | 6 |
| Terumo Hypodermic Syringes with Needle | 6 |
| Terumo Hypodermic Syringes without Needle | 6 |
| Terumo Pharmacy Trays | 7 |

Blood Collection Sets

| | |
|---|---|
| Surshield® Safety Blood Collection Sets | 8 |
| Multi-Sample Luer Adapter | 8 |

Winged Infusion Needle Sets

| | |
|--|---|
| Surshield® Safety Winged Infusion Sets | 9 |
| SURFLO® Winged Infusion Sets | 9 |

IV Catheters

| | |
|---------------------------------------|----|
| SURFLO® ETFE I.V. Catheters | 10 |
| SurFlash® Polyurethane I.V. Catheters | 10 |

Frequently Asked Questions and Technical Support

11



Pick how you click.




With 3 ways to activate and a broader range of sizes than most other brands, SurGuard®3 is the industry's most adaptable safety needle. Add separate thumb and finger grips, sharper needles, and competitive pricing, and you'll understand why we say SurGuard®3 is the one needle that fits *every* clinician's style.

- 3 modes of activation—thumb, finger, or hard surface
- Separate finger and thumb grips with ridges and nonslip surfaces
- Hub and needle locks for security post activation
- Broad range of sizes
- Firm injection feel
- Sharpest needles among major brands tested*
 - 10% sharper than similar needles from Becton, Dickinson and Company
 - More comfortable injections for patients

Bottom line: No matter what you're looking for in a needle, you're looking for SurGuard®3.

*Reference: Data on file. Terumo Medical Products, April 2009.

Order the industry's most adaptable safety needle today!



SurGuard®3 Safety Hypodermic Needle

Packaging: 100/box, 8 boxes/case (800)

| Order No. | Description |
|-----------|---------------------|
| SG3-1825 | 18G x 1" needle |
| SG3-1838 | 18G x 1 1/2" needle |
| SG3-2025 | 20G x 1" needle |
| SG3-2038 | 20G x 1 1/2" needle |
| SG3-2051 | 20G x 2" needle |
| SG3-2125 | 21G x 1" needle |
| SG3-2138 | 21G x 1 1/2" needle |
| SG3-2225 | 22G x 1" needle |
| SG3-2238 | 22G x 1 1/2" needle |
| SG3-2325 | 23G x 1" needle |
| SG3-2338 | 23G x 1 1/2" needle |
| SG3-2516 | 25G x 5/8" needle |
| SG3-2525 | 25G x 1" needle |
| SG3-2538 | 25G x 1 1/2" needle |
| SG3-2713 | 27G x 1/2" needle |
| SG3-3013 | 30G x 1/2" needle |

SurGuard®3 Hypodermic Syringes with Safety Needle

Packaging: 100/box, 4 boxes/case (400)

| Order No. | Description |
|-------------|--------------------------------------|
| SG3-01T2516 | 1CC Syringe with 25G x 5/8" needle |
| SG3-01T2713 | 1CC Syringe with 27G x 1/2" needle |
| SG3-03L2025 | 3CC Syringe with 20G x 1" needle |
| SG3-03L2038 | 3CC Syringe with 20G x 1 1/2" needle |
| SG3-03L2125 | 3CC Syringe with 21G x 1" needle |
| SG3-03L2138 | 3CC Syringe with 21G x 1 1/2" needle |
| SG3-03L2225 | 3CC Syringe with 22G x 1" needle |
| SG3-03L2238 | 3CC Syringe with 22G x 1 1/2" needle |
| SG3-03L2325 | 3CC Syringe with 23G x 1" needle |
| SG3-03L2516 | 3CC Syringe with 25G x 5/8" needle |
| SG3-03L2525 | 3CC Syringe with 25G x 1" needle |



You know what you like.

And *whatever* your preference, you'll find it in SurGuard®3!

SurGuard®3
SAFETY HYPODERMIC NEEDLE

Our double-bevel design is precision-ground and honed to exacting quality standards. This attention to detail enables Terumo needles to penetrate skin with less force, which can mean greater comfort for patients.




Terumo® Standard Hypodermic Needles

Packed: 100/box, 10 boxes/case (1,000)

| Product Code | Description |
|--------------|-------------------|
| ■ NN1825R | 18G x 1" needle |
| ■ NN1838R | 18G x 1½" needle |
| ■ NN1925R | 19G x 1" needle |
| ■ NN1938R | 19G x 1½" needle |
| ■ NN2025R | 20G x 1" needle |
| ■ NN2038R | 20G x 1½" needle |
| ■ NN2125R | 21G x 1" needle |
| ■ NN2138R | 21G x 1½" needle |
| ■ NN2225R | 22G x 1" needle |
| ■ NN2238R | 22G x 1½" needle |
| ■ NN2325R | 23G x 1" needle |
| ■ NN2338R | 23G x 1½" needle |
| ■ NN2516R | 25G x 5/8" needle |
| ■ NN2525R | 25G x 1" needle |
| ■ NN2538R | 25G x 1½" needle |
| ■ NN2613R | 26G x ½" needle |
| ■ NN2713R | 27G x ½" needle |
| ■ NN2732R | 27G x 1¼" needle |
| ■ NN3013R | 30G x ½" needle |
| ■ NN3025R | 30G x 1" needle |

All Competitors



* In Independent laboratory testing, Terumo needles were shown, on average, to be the sharpest among brands tested:

10% sharper than BD (22G - 30G)

16% sharper than Covidien (22G - 27G)

These products have no components made of natural rubber latex.

*Data on file. Terumo Medical Products, April 2016.

Conventional Syringes

1cc Tuberculin Syringes

Graduations: 0.01cc (1cc total)

TERUMO TB Syringes have clear, bold scale markings which assure dosage control and accuracy. The package is tamper evident, for greater protection of the sterile contents.

Packed: 100/box, 10 boxes/case (1,000)

Product Code Description

| | |
|------------|---|
| SS-01T | 1cc TB Syringe without needle, with luer slip |
| SS-01T2516 | 1cc TB Syringe with 25G x 5/8" removable needle |
| SS-01T2609 | 1cc TB Syringe with 26G x 3/8" removable needle |
| SS-01T2713 | 1cc TB Syringe with 27G x 1/2" removable needle |



Catheter Tip



Luer Lock Tip



Eccentric Luer Tip



Luer Slip Tip

3cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

| Product Code | Description |
|--------------|--------------------------|
| SS-03L2025 | 3cc, 20G x 1" needle |
| SS-03L2038 | 3cc, 20G x 1 1/2" needle |
| SS-03L2125 | 3cc, 21G x 1" needle |
| SS-03L2138 | 3cc, 21G x 1 1/2" needle |
| SS-03L2225 | 3cc, 22G x 1" needle |
| SS-03L2238 | 3cc, 22G x 1 1/2" needle |
| SS-03L2325 | 3cc, 23G x 1" needle |
| SS-03L2338 | 3cc, 23G x 1 1/2" needle |
| SS-03L2516 | 3cc, 25G x 5/8" needle |
| SS-03L2525 | 3cc, 25G x 1" needle |

5cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.2cc (6cc total)

Packed: 100/box, 6 boxes/case (600)

| Product Code | Description |
|--------------|--------------------------|
| SS-05L2025 | 5cc, 20G x 1" needle |
| SS-05L2038 | 5cc, 20G x 1 1/2" needle |
| SS-05L2125 | 5cc, 21G x 1" needle |
| SS-05L2138 | 5cc, 21G x 1 1/2" needle |
| SS-05L2225 | 5cc, 22G x 1" needle |
| SS-05L2238 | 5cc, 22G x 1 1/2" needle |

10cc Syringes with Needle

Style: Luer Lock Tip

Graduations: 0.2cc (12cc total)

Packed: 100/box, 6 boxes/case (600)

| Product Code | Description |
|--------------|---------------------------|
| SS-10L2025 | 10cc, 20G x 1" needle |
| SS-10L2038 | 10cc, 20G x 1 1/2" needle |
| SS-10L2125 | 10cc, 21G x 1" needle |
| SS-10L2138 | 10cc, 21G x 1 1/2" needle |

3cc Syringes without Needle

Graduations: 0.1cc (3cc total)

Packed: 100/box, 10 boxes/case (1,000)

| Product Code | Description |
|--------------|---------------|
| SS-03L | Luer Lock Tip |
| SS-03S | Luer Slip Tip |

5cc Syringes without Needle

Graduations: 0.2cc (6cc total)

Packed: 100/box, 6 boxes/case (600)

| Product Code | Description |
|--------------|-------------------|
| SS-05L | 5cc Luer Lock Tip |
| SS-05S | 5cc Luer Slip Tip |

10cc Syringes without Needle

Graduations: 0.2cc (12cc total)

Packed: 100/box, 6 boxes/case (600)

| Product Code | Description |
|--------------|-------------------------|
| SS-10L | 10cc Luer Lock Tip |
| SS-10S | 10cc Slip Tip |
| SS-10ES | 10cc Eccentric Luer Tip |

20cc Syringes without Needle

Graduations: 1.0cc (25cc total)

Packed: 50/box, 4 boxes/case (200)

| Product Code | Description |
|--------------|-------------------------|
| SS-20L2 | 20cc Luer Lock Tip |
| SS-20S2 | 20cc Luer Slip Tip |
| SS-20ES2 | 20cc Eccentric Luer Tip |

30cc Syringes without Needle

Graduations: 1.0cc (35cc total)

Packed: 25/box, 8 boxes/case (200)

| Product Code | Description |
|--------------|-------------------------|
| SS-30L | 30cc Luer Lock Tip |
| SS-30S | 30cc Luer Slip Tip |
| SS-30ES | 30cc Eccentric Luer Tip |

60cc Syringes without Needle

Graduations: 1.0cc (60cc total)

Packed: 25/box, 4 boxes/case (100)

| Product Code | Description |
|--------------|-------------------------|
| SS-60L | 60cc Luer Lock Tip |
| SS-60ES | 60cc Eccentric Luer Tip |
| SS-60C | 60cc 2oz. Catheter Tip |

These products have no components made of natural rubber latex.

Pharmacy Trays

Terumo® Pharmacy Trays

- Terumo syringes are neatly stacked inside a tray for easy removal.
- Convenient peel-back cover insures sterility.
- Trays are individually embossed with the appropriate syringe size for accurate identification when stacked on shelves.
- The unique tray design requires no additional packaging components, which reduces waste and increases cost efficiency.



| Product Code | Description |
|--------------|--|
| SS10LT | 10cc Luer Lock Tip without Needle Pharmacy Tray, Graduations: 0.2cc (12cc total) Packed: 25 syringes/tray, 8 trays/case (200) |
| SS20LT | 20cc Luer Lock Tip without Needle Pharmacy Tray, Graduations: 1.0cc (25cc total) Packed: 25 syringes/tray, 8 trays/case (200) |
| SS30LT | 30cc Luer Lock Tip without Needle Pharmacy Tray, Graduations: 1.0cc (35cc total) Packed: 10 syringes/tray, 12 trays/case (120) |
| SS60LT | 60cc Luer Lock Tip without Needle Pharmacy Tray, Graduations: 1.0cc (60cc total) Packed: 10 syringes/tray, 12 trays/case (120) |

These products have no components made of natural rubber latex.



Surshield® Safety Winged Blood Collection Sets

Packed: 50/box, 10 boxes/case (500)

Product Code **Description**

(12" Tubing)

| | |
|-------------|-------------------|
| MN*SVS21B30 | 21G x 3/4" needle |
| MN*SVS23B30 | 23G x 3/4" needle |
| MN*SVS25B30 | 25G x 3/4" needle |

(7" Tubing)

| | |
|-------------|-------------------|
| MN*SVS21B18 | 21G x 3/4" needle |
| MN*SVS23B18 | 23G x 3/4" needle |



Multi-Sample Luer Adapter

Sterile and disposable, adapters fit into standard infusion sets transforming them into winged blood collection sets. This allows multiple draws into evacuated tubes. Accepts standard tube holders.

Packed: 100/box, 30 boxes/case (3,000)

Product Code **Description**

| | |
|-------------|---|
| XX*MN2000TE | VENOJECT® Multi-Sample Luer Adapter (male). |
|-------------|---|




Blood Collection Sets

Safety


Terumo's Surshield Safety Winged Needle Sets with one-handed activation provides a solution with a cost-effective, easy to use design. Press shield over the needle until it locks into place.

- Unique needle locking mechanism helps reduce the risk of needlestick injuries
- One-handed activation for safety and ease of use
- No change in needle insertion technique
- Simple, built-in safety design is easy to use

1. Flip safety shield toward needle. Remove the needle and apply digital pressure to the site.



2. With the wing and shield between your thumb and index finger, press the safety shield against a hard surface such as a bedside table until an audible click is heard.



3. Visually confirm activation of the safety feature. (10/08)



Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

These products have no components made of natural rubber latex.



Surshield® Safety Winged Infusion Sets

Packed: 50/box, 10 boxes/case (500)

| Product Code | Description |
|---------------------|-----------------|
| (3½" Tubing) | |
| ■ SV*S19BLS | 19G x ¾" needle |
| ■ SV*S21BLS | 21G x ¾" needle |
| ■ SV*S23BLS | 23G x ¾" needle |
| ■ SV*S25BLS | 25G x ¾" needle |

| | |
|---------------------|-----------------|
| (12" Tubing) | |
| ■ SV*S19BL | 19G x ¾" needle |
| ■ SV*S21BL | 21G x ¾" needle |
| ■ SV*S23BL | 23G x ¾" needle |
| ■ SV*S25BL | 25G x ¾" needle |



SURFLO® Winged Infusion Sets

| Product Code | Description |
|---------------------|-------------|
| (3½" Tubing) | |

Packed: 100/box, 5 boxes/case (500)

| | |
|------------|-----------------|
| ■ SV*19BLS | 19G x ¾" needle |
| ■ SV*21BLS | 21G x ¾" needle |
| ■ SV*23BLS | 23G x ¾" needle |
| ■ SV*25BLS | 25G x ¾" needle |

| | |
|---|--|
| (8" Tubing) | |
| Packed: 100/box, 10 boxes/case (1,000) | |

| | |
|-----------|-----------------|
| ■ SV*25EL | 25G x ½" needle |
| ■ SV*27EL | 27G x ½" needle |

| | |
|--|--|
| (12" Tubing) | |
| Packed: 50/box, 10 boxes/case (500) | |

| | |
|------------|-----------------|
| ■ SV-18BLK | 18G x ¾" needle |
| ■ SV-19BLK | 19G x ¾" needle |
| ■ SV-21BLK | 21G x ¾" needle |
| ■ SV-22BLK | 22G x ¾" needle |
| ■ SV-23BLK | 23G x ¾" needle |
| ■ SV-25BLK | 25G x ¾" needle |

I.V. Catheters



Terumo's Polyurethane I.V. Catheter enables quick and easy confirmation of placement in the vein.

With SurFlash, visual confirmation of catheter placement occurs due to a unique flashback feature. Flashback appears between the specially grooved needle and the catheter, providing medical staff with a higher certainty of proper catheter placement. The specially grooved needle also improves the ability to visualize flashback in patients where venipuncture is more difficult, and because SurFlash is made of polyurethane, which is firm at room temperature and softens at body temperature, patient comfort is enhanced, without compromising proper technique.

SURFLO® ETFE I.V. Catheters (Large Photo)

A sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue. Preview chamber in clear sure-grip hub gives immediate indication of proper placement. Flexible thinwall catheter with large inner diameter assures good blood flow. All this plus safe, easy-to-handle, tamper-evident packaging make SURFLO® I.V. Catheters a great clinical value.

Packed: 50/box, 4 boxes/case (200)

| Product Code | Description |
|--------------|-------------|
| SR-OX1451CA | 14G x 2" |
| SR-OX1632CA | 16G x 1¼" |
| SR-OX1651CA | 16G x 2" |
| SR-OX1664CA | 16G x 2½" |
| SR-OX1832CA | 18G x 1¼" |
| SR-OX1851CA | 18G x 2" |
| SR-OX1864CA | 18G x 2½" |
| SR-OX2025CA | 20G x 1" |
| SR-OX2032CA | 20G x 1¼" |
| SR-OX2051CA | 20G x 2" |
| SR-OX2225CA | 22G x 1" |
| SR-OX2419CA | 24G x ¾" |

SurFlash® Polyurethane I.V. Catheters (Not Pictured)

- Visual confirmation of catheter placement during venipuncture occurs in a flash!
- Made of polyurethane, which makes it extremely flexible, kink-resistant and durable
- Provides medical staff with higher certainty of catheter placement due to unique flashback feature
- Improves visualization for difficult venipunctures

Packed: 50/box, 4 boxes/case (200)

| Product Code | Description |
|--------------|-------------|
| SR*FF1451 | 14G x 2" |
| SR*FF1651 | 16G x 2" |
| SR*FF1664 | 16G x 2 ½" |
| SR*FF1832 | 18G x 1 ¼" |
| SR*FF1851 | 18G x 2" |
| SR*FF1864 | 18G x 2 ½" |
| SR*FF2025 | 20G x 1" |
| SR*FF2032 | 20G x 1 ¼" |
| SR*FF2051 | 20G x 2" |
| SR*FF2225 | 22G x 1" |
| SR*FF2419 | 24G x ¾" |

These products have no components made of natural rubber latex.

Frequently Asked Questions and Technical Support

Who is the Sales Rep for my area?

Please contact us by phone or email with your address to obtain the name of your local Territory Manager. You can also visit our website, www.TerumoTMP.com, for more information.

How can I obtain samples?

Please contact your local Territory Manager or visit our website to request samples.

How can I order Terumo products through a distributor?

Please contact your authorized distributor. To find a distributor, visit us at www.TerumoTMP.com/Support to find a List of Authorized Terumo Distributors. Distributors should call Customer Administration for all questions related to orders and invoices: 800.888.3786

How can I obtain pricing?

You must check with your local distributor for pricing since our business is conducted through authorized distributors.

Where are you located?

Terumo Medical Canada, Inc.
10911 Keele Street
Units 2-4
Vaughan, Ontario
L6A 5A6
Phone: 833.883.7866
Fax: 647.715.5973
terumocanadacustomer.admin@terumomedical.com

What is the shelf life of Terumo Products?

Standard & Safety Syringes & Needles: 5 years
Standard IV Catheters: 3 years
Standard & Safety Blood Collection/Infusion: 3 years

What's the difference between SurFlash® & Surflo® IV Catheters?

SurFlash is an I.V. Polyurethane Catheter which allows for a visual confirmation of catheter placement due to its unique flashback feature. SURFLO is our regular ETFE I.V. Catheter. Both are radio-opaque and these products have no components made of natural rubber latex.

www.TerumoTMP.com

RX ONLY. The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURGUARD is a registered trademarks of Terumo Medical Corporation, Somerset, NJ. SURFLO, SURSHIELD and VENOJECT are registered trademarks of Terumo Corporation, Tokyo, Japan and they are registered with the U.S. Patent and Trademark Office. Approval # PM-03466