

# Sur-Vet® Surflo® I.V. Catheters

A sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue.

Preview chamber in clear sure-grip hub gives immediate indication of proper placement.

Flexible thinwall catheter with large inner diameter assures good blood flow.

All this plus easy-to-handle, tamper-evident packaging make Sur-Vet® SURFLO® I.V. Catheters a great clinical value.



## Sur-Vet® Surflo ETFE I.V. Catheters

Packed: 50/box, 4 boxes/case (200)

Product Code	Description	Needle Gauge
SR-OX1832V	18G x 1 ¼"	20G
SR-OX2025V	20G x 1"	22G
SR-OX2032V	20G x 1 ¼"	22G
SR-OX2051V	20G x 2"	22G
SR-OX2225V	22G x 1"	24G
SR-OX2419V	24G x ¾"	27G

## SURFLO® ETFE I.V. Catheters

Packed: 50/box, 4 boxes/case (200)

Product Code	Description	Needle Gauge
SR-OX1451CA	14G x 2"	16G
SR-OX1632CA	16G x 1 ¼"	18G
SR-OX1651CA	16G x 2"	18G
SR-OX1664CA	16G x 2 ½"	18G
SR-OX1851CA	18G x 2"	20G
SR-OX1864CA	18G x 2 ½"	20G



Call your Terumo representative at **1.800.888.3786** or visit [www.TerumoTMP.com](http://www.TerumoTMP.com) and click on Find Your Local Terumo Representative.



# Bred for Performance

These products contain no components made of natural rubber latex.

**RX ONLY.** The advertisement is directed to veterinarians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SUR-VET and SURFLO are registered trademarks of Terumo Corporation, Tokyo, Japan and they are registered with the U.S. Patent and Trademark Office. © 2023 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. Rev. 1/23 PM-06589